KONUSCIENCE

100X-1200X ZOOM MICROSCOPE

www.konus.com

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
BATTERY INSTALLATION

- 1) Remove the LED lamp from the arc bracket.

- 2) With a small Phillips screwdriver, loosen the screw on the side.

- 3) Carefully remove the mirror cover, make sure the glass mirror does not detach from the plastic frame.

- 4) Discard the old batteries.

- 5) Install two 1.5 Volt LR41/AG3 batteries per the polarity markings in the battery compartment.


- 6) Replace the mirror cover and tighten the screw.

- 7) Install the LED lamp back to the bracket.

PROCEDURES

- 1) Tilt the body and adjust the location of the reflector. Therefore, the light is fully caught by the mirror.

- 2) If the experiment room is dim with low intensity of light, or the focus is unclear at extreme magnification, it is recommended to turn on the reflector. Automatically the light source lights will allow for observation.

- 3) When the light is fully reflected by the mirror, which can be seen via the eyepiece, the microscope will be ready for inspection.

- 4) Now, put the prepared slide on the stage, and fasten it in place with the clips.

- 5) Next, close which magnification strength you desire. Remember the greater the length of the objective lens, the bigger the magnification. In general, inspection is usually made at low setting. To change the magnification strength

turn the revolving turret until you hear a click.

- 7) Using the focusing knob, let down the lens as close as possible to the prepared slide without making interactions. Next, looking through the eyepiece, turn the knob anticlockwise until the reflection achieves clarity.

HOW TO CHANGE THE ILLUMINATION OF THE MAGNIFICATION FIELD

If you are working with the lowest power objective (reduced magnification) the light on the field will probably appear too strong and dazzling. You can solve your problem by rotating the disk diaphragm located beneath the stage: through one of its smallest apertures, you will receive the appropriate illumination. Very thin, low-coloured and almost transparent specimens will consequently exploit high-contrast details coming from the use of a smaller diaphragm aperture allowing easier observations.


In case the field of view is too dark, such as when you are using high power objectives, first of all make sure that the illuminator did not move from its original position and then check that you are using the wider apertures of the disk diaphragm. Sometimes the specimen can look dark because of its thickness, do not forget that with the optical microscope you can examine only thin and transparent specimen through which light can pass, if not, it will be impossible to see the smallest details.

EYEPIECES, OBJECTIVES AND MAGNIFICATIONS

Eyepieces and objectives are lens systems enlarging images. The magnification at which you are observing the specimen is simply to calculate, as it is the result between the eyepiece magnification multiplied for the magnification of the objectives. Objectives are much more perfect systems, because they have the hardest magnification task. The advantage coming from a rotatable nosepiece holding three objectives, is that you can easily change the magnifications by simply rotating the nosepiece. The standard eyepiece in your microscope offers 10 or 20 magnifications, as it is a zoom eyepiece, it is necessary to rotate the knurled ring to change the zoom power of the eyepiece from 10 to 20.


HOW TO USE THE CAMERA ADAPTER

The camera adapter is used to apply the camera on the microscope and take pictures using the microscope. It is composed of two parts. Place the base, on the microscope eyepiece tube, insert the camera paid attention that its objective will be exactly on the microscope eyepiece, now, fix the camera using the plastic adjustments.

USE GUIDANCE OF VIEWER

- 1) Turn left the eyepiece for microscope to dislocate the eyepiece and insert the viewer head in the microscope proper. In the case of zoom microscope, turn left the nut underneath the eyepiece and insert the viewer head in the body of microscope and turn right the nut tightly.


Put the plastic viewer screen into the hole given in the right (or upper) front side of viewer head. The viewer screen is set with its coarse face positioned inside.

- 2) Set prepared slide in the clips of stage. Make adjustment so that the light from the lamp may be perfectly cast the stage hole onto the objective lens.

Then operate the handle to move the tube up and down, which causes the prepared slide to be projected on the viewer screen.

When the image produced on the viewer screen is not clear, adjust the illumination direction of the light source lamp. (Do not turn over 360° the section the light source lamp is fixed on the arm. Otherwise it will cause a breakdown).


- 3) When observation is made with viewer, the magnification is about 1.2 times that of

objective lens. Observation by the use of viewer can be made without interruption for about one hour because of the life time of battery. Please save battery life by putting off it while not in use.


HOW TO USE THE PROJECTION DEVICE

- 1) Turn the eyepiece to the left to dislocate it, and insert the viewer head in the microscope properly.
- 2) Face it toward white paper or a white wall that you are able to see projected image. Otherwise, attach a piece of white paper on the wall if the wall colour is dim.
- 3) Unlit the room. Mount a prepared slide on the stage.
- 4) Affix the condenser lens cap to the illumination lamp. Then turn the light on.
- 5) Revolve the viewing projection and drawing head towards the wall at a distance of 1.5m from the wall.
- 6) Adjust the illumination lamp so that it is as bright as possible.
- 7) Focus the image by carefully erecting the objective.
- 8) Once in focus, adjust the illumination lamp one more time.
- 9) You can see a clear image projected. To increase the magnification, move the microscope away from the wall.


viewing screen, darken the room by turning off the light.

- 3) Place a piece of white paper horizontally in front of the foot of the microscope.
- 4) Station the reflecting mirror and rotate the focusing knob until a pleasing image is projected.


HOW TO USE THE HAND MICROTOME

With a little effort and a lot of satisfaction, you can prepare on your own very thin sections of vegetables. Take your microtome (fig.11), a tool you will find among the accessories of your microscope. Make sure the diameter of your material can fit one of the two comma-shaped holes of the microtome. Insert one end of your vegetable in one of the two holes and let a little portion stick out. Keeping the vegetal material pressed against the narrowest part of the comma-shaped hole, rotate the blade towards the vegetable, until, under a uniform pressure, you cut it. The section must be as thin as possible: the thinner, the lesser the vegetable portion will be sticking out of the comma-shaped hole of the microtome.

HOW TO MAKE USE OF THE DRAWING DEVICE

- 1) Position and adjust the arm to a vertical position.
- 2) After the image has been projected on the

HOW TO MAKE A PREPARED SLIDE: • TEMPORARY MOUNT

- Wipe the slide and cover glass clean
- Thin the sample with the microtome (be extremely careful)
- Then pick up with tweezers, and put it on the central part of the glass slide.
- Next, add one drop of water on the sample with a dissecting needle. If the sample is clear, add one drop of methylene blue or eosin solution (be extremely careful).
- Gently put the cover glass on it, take care not to let any air bubbles in it.
- Remove any excess water or dyeing solutions with blotting paper.
- Now, it is ready for observation.

• PERMANENT MOUNT

- Wipe the slide and cover glass clean
- Continue as above, but before covering the slide with the cover glass, add few drops of gum media (or Canada Balsam) solution or transparent adhesive glue with a dissecting needle to the slide.
- Put down on the cover glass with tweezers or fix it in place, and leave it to dry for about a day.

SPECIMENS READY TO USE

At the same shop where you bought your microscope, Konus make available ready-to-use specimen series of different kind. Such didactical series consists of material which you would never be able to prepare on your own, either because it would be impossible to obtain it (who indeed would be able to give you a piece of human liver?), or because it would be too difficult to prepare it with the same expertness of specialized scientific laboratories.

Konus collection consists of didactical series helpful to consider very different, updated and interesting topics. The line includes several series, each one of them packed and presented in an elegant case also containing an introductory note explaining the series features.

THE COMPLETE PROGRAM OF KONUS MICROSCOPE SLIDES SERIES

Each series contains 10, 15 or 25 microscope slides of best quality and an explanatory text:

- #4960 Series for the teaching of science
- #4963 Series for the teaching of biology (Part I)
- #4964 Series for the teaching of biology (Part II)
- #4976 General Biology
- #4977 Zoology, invertebrate and insects
- #4978 Zoology, vertebrate including mammalia
- #4886 The wonderful structure of vegetable cells
- #4887 Cell development and reproduction

- #4881 Plant reproduction and propagation
- #4900 The wonderful world in a drop of water
- #4950 The life of the soil
- #4986 Biology, the cell and animal tissues
- #4910 The human body: normal tissues
- #4911 The human body: normal tissues
- #4987 Human body, the structure of organs
- #4914 The human body: pathological tissues -1
- #4915 The human body: pathological tissues 2
- #4855 Fish, Frogs and other Amphibians
- #4856 Lizards, Snakes and Birds
- #4979 Zoology, Protozoa
- #4980 Zoology, Insects (25 slides)
- #4981 Zoology, Reproduction in invertebrates
- #4862 Animal reproduction
- #4863 Animal embryology and development
- #4918 Human and animal parasites
- #4877 Bacteria and simple organisms
- #4982 Botany, Fungi, Lichens and Briophyta
- #4983 Botany, Hepaticae, Mosses, Pteridophytae & Gymnosperms
- #4984 Botany, angiospermae, monocotyledons (Liliopsidae)
- #4985 Botany, angiospermae, dicotyledons (magnoliopsidae)
- #4952 Animals and plants damaged by the environment

HOW TO KEEP YOUR MICROSCOPE IN GOOD CONDITIONS

Before you put the microscope away, make sure you have removed the slide from the stage. Do not forget to clean it and take off the dust with a brush using a soft cloth for the eyepiece. Eyelashes in fact frequently deposit sebum on the upper lens of the eyepiece, thus worsening your observations. After you have vigorously blown on the eyepiece breath on the upper lens, and, through a soft cloth, clean the lens following a rotating direction. Remember that if you insist too strongly on it, especially if it covered with dust, it will become like a frosted glass. Objectives do not get dirty as much as the eyepiece does. Should they need to be cleaned, take a wad of cotton-wool would on a toothpick and gently rotate it on the lens. Never use detergents or, even worse alcohol and similar products, because you may melt the plastic parts of the lens mount or the paint may come off irretrievably damaging the lenses. At last, remember to remove batteries from the base of the microscope, so that they will not run down when you do not use the microscope. It may occur, for instance that the illuminator, switches the contacts. Batteries can run for about one hour.

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© KONUS 2020 ALL RIGHT RESERVED - PRINTED IN P.R.C.
KONUS IS REGISTERED TRADE MARK OF KONUS IN U.S.A., ITALY, SPAIN AND SELECTED