|     |                                  |    |
|-----|----------------------------------|----|
| D   | Elektrische Zahnbürste           |    |
| NL  | Elektrische tandenborstel        |    |
| F   | Brosse à dents électrique        |    |
| E   | Cepillo de dientes eléctrico     |    |
| I   | Spazzolino da denti elettrico    |    |
| GB  | Electric Toothbrush              |    |
| PL  | Elektryczna szczoteczka do zębów |    |
| H   | Elektronos fogkefe               |    |
| UA  | Електрична зубна щітка           |    |
| RUS | Электрическая зубная щётка       |    |
| AR  | فريسة الأسنان الكهربائية         |    |
| D   | Bedienungsanleitung/Garantie     | 02 |
| NL  | Gebruiksaanwijzing               | 03 |
| F   | Mode d'emploi                    | 04 |
| E   | Instrucciones de servicio        | 05 |
| I   | Istruzioni per l'uso             | 06 |
| GB  | Instruction Manual               | 07 |
| PL  | Instrukcja obsługi/Gwarancja     | 08 |
| H   | Használati utasítás              | 09 |
| UA  | Інструкція з експлуатації        | 10 |
| RUS | Руководство по эксплуатации      | 11 |
| AR  | فرشاة الأسنان الكهربائية         | 12 |

## ELEKTRISCHE ZAHNBÜRSTE

EZ 5622

**CE**

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### 2 Deutsch

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## 7 English

### User manual

Thank you for choosing our product. We hope you will enjoy using the appliance.

### Safety instructions

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- Do not immerse the device in water or other liquids.
- Do only use original parts.
- This device serves for the dental and oral hygiene of adult persons.
- Children aged 7 and above may use the device under the supervision of an adult.

**WARNING:**  
This device is not suited for children under the age of 3 as small parts could be broken off and swallowed.

- Store the device out of the reach of children. Children may not play with the device.

### Handling of batteries

**⚠ WARNING: Danger of explosion!**  
Do not expose the batteries to high temperatures or direct sunlight. Never throw batteries into fire.

**⚠ CAUTION:**  
Different battery types or new and spent batteries may not be used together.

### Inserting the batteries

(Battery is not included in the scope of delivery)

- Open the battery compartment. Pull off the cover of the stand.
- Insert 2 batteries of the type AA/R6. Ensure the correct polarity! You will find a marking on the contacts inside the device. Set the positive pole of the battery onto the positive contact in the device and the negative pole of the battery onto the negative contact in the device.
- Close the battery compartment again.

**NOTE:**

- If the device provides too little power or if it shuts off during use, replace the batteries against new ones.
- You may also insert charged rechargeable batteries of the AA type.
- If the device is not used for a longer period of time, remove the batteries/ rechargeable batteries to prevent a "leaking out" of battery acid.

### Usage

- Stick the brush attachment onto the device. The brush head will show in the direction of the **I / O** switch. Remove the protective cap from the brush head, if necessary.
- Moisten the brush head and apply some toothpaste.
- Rinse your mouth with water.
- Guide the brush to your teeth before turning on the device with the I switch. You will thus prevent splatters.
- Place the brush against the edge of the gums. Brush the outsides, the insides and the occlusal surface of your teeth. Slowly guide the brush head from tooth to tooth. Only exert little pressure. Dentists recommend a brushing period of a total of 2 minutes.
- Shut off the device with the **O** switch.

**NOTE:**

- During the first few uses it can result in a slight bleeding of the gums.
- If these symptoms should last for longer than 2 weeks, you should contact a dentist.

### Cleaning

**CAUTION:**

- Do not immerse the device in water. This could destroy the electronics.
- Do not use sharp or abrasive cleaning agents.

- Switch off the appliance before cleaning.
- You may wipe the housing of the device with a slightly damp cloth.
- Pull off the brush attachment. Clean it under running water after each use.

### Technical Specifications

Model: .....EZ 5622  
Net weight: .....0,80 kg  
Battery: .....2x 1,5 V Type "AA" / "R6"

The right to make technical and design modifications in the course of continuous product development remains reserved.

This device complies with all current CE directives, such as electromagnetic compatibility and low voltage directive and is manufactured according to the latest safety regulations.

### Disposal

- Remove the battery from the battery compartment.
- Dispose of the battery at an appropriate collection point or distributor. This product must not be disposed of in standard domestic refuse.
- Do not dispose of the device with household waste at the end of its service life. Bring it to an official collection point for recycling. By doing so, you help protect the environment.

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