ORDO

HYDRO SONIC

WATER FLOSSER

Important Information

Wichtige Informationen

Information Capitale

Información importante

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Viktig information

Vigtig information

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Oluline teave

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ordolife.com

Ordo, Manchester, WA14 4DR, UK.

Authorised Rep Compliance Ltd., Ground Floor,
71 Lower Baggot Street, Dublin, D02 P593, Ireland.

English

Please read all usage instructions carefully before first use and save for future reference.

Technical Data

Model number: HSWF01-CG, HSWF01-RG, HSWF01-W, HSWF01-PV, HSWF01-IB

Input: DC 5V, 2.5W

Protection class: III

IP rating: IPX7

Dimensions, with no floss tip installed (WxDxH): 7.4 x 9.5 x 23.1 cm

Dimensions, with Premium Clean floss tip installed (WxDxH): 7.4 x 9.5 x 30.6 cm

Rechargeable battery: lithium-ion, 3.7V, 1500mAh

Charging time: 6 hours

Charging and operating temperature between: 5°C and 40°C

Battery life: 30 days/full charge (at 2 minutes of use per day).

Water tank capacity: 260ml, providing 55 seconds of use.

Timer: deactivation after 2 minutes of operation.

Warranty

This product is covered by a 2-year warranty. The extent of the warranty on this product is limited to technical defects caused by faulty production processes. In the event you want to claim warranty, please be sure to contact our customer service team for instructions. For more information, please visit ordolife.com/warranty

60-Day Money-Back Guarantee

Satisfaction guaranteed – if you don't love your water flosser within 60 days, we'll give you a refund. For more information on our guarantee please visit ordolife.com/guarantee

Customer Service

For more information on Ordo's products and services please visit ordolife.com

If the device is damaged, defective, in need of repair or you require assistance please contact our customer service team:

Email: hello@ordolife.com

Support: ordolife.com/help

Our team are always happy to help.

Recycling & Disposal

At Ordo we are always looking for solutions to reduce packaging and product materials that cannot be recycled effectively. We ask every member of the Ordo family to help us by recycling responsibly.

For more information, please visit ordolife.com/recycling

Do not dispose of the product with normal household waste at the end of its life. Take it to an official collection point for recycling and proper disposal.

This product contains a built-in, lithium-ion rechargeable battery. The battery must only be removed by a qualified professional when the appliance is disposed of.

WARNINGS & CAUTIONS

Intended use

The water flosser is designed exclusively for cleaning teeth and gums. It is only intended for private use and is not suitable for commercial purposes. Only use the water flosser as described in the operating instructions. Any other use is deemed improper and may result in personal injury or damage to property. The water flosser is not a child's toy. Ordoclub Ltd accepts no liability for injury or damage caused by improper or incorrect use.

General warnings

- Only use original Ordo accessories and consumables. This includes the floss tips and USB cable.
- This device cannot be serviced or repaired by the user. If the device is damaged, stop using it and contact our customer service team.
- Liability and warranty claims are waived in the event of repairs performed by the user, improper connection of the device or incorrect operation.

Risk of electric shock

A faulty electrical installation or excessive mains voltage may result in an electric shock.

- Only connect the device to a power source if its output voltage corresponds to the specification on the device.
- Only connect the device to an easily accessible socket so that you can quickly disconnect it from the power source in case of a failure.
- Do not operate the device if it exhibits visible damage, e.g. cracks, breaks or deformations.

Check the device for damage from time to time.

- Do not use the USB cable if it exhibits visible damage or if the main cord or USB connector is defective. Check the USB cable for damage from time to time. The USB cable cannot be repaired if damaged. In the event of damage, the USB cable must be disposed of and replaced.
- Do not disassemble or modify the device.
- The device must not be charged with an external timer or with a separate tele-control system.
- Do not immerse the USB cable in water or other liquids.
- Never touch the USB connectors or USB-C port with damp hands.
- Ensure the entire device and the USB-C input are fully dry before charging the device.
- Only charge the device in rooms located indoors.
- Never charge the device such that there is a risk of it falling into a bathtub or sink.
- Never attempt to grab the device if it has fallen into water during charging. In such a case, immediately disconnect the USB cable, if it is safe to do so.
- Only charge the device with an original Ordo USB cable, one of which is enclosed with the product.
- Disconnect the USB cable when you are not charging the device, before you clean the device or in the event of a fault.
- Never unplug the USB cable by pulling on the main cord, use the plastic casing of the USB connector instead.
- Ensure the USB-C input cover on the device is securely closed when the device is not being charged.
- Never use the USB cable as a carrying strap.
- Do not kink the USB cable and do not lay it over sharp edges.
- Do not insert foreign objects into any of the USB connectors or any opening on the device.

Battery safety information

Improper handling of rechargeable batteries may cause an explosion. Fluid leaked from the rechargeable battery may lead to chemical burns.

- Only use the rechargeable battery provided with the device.
- Never short circuit, overcharge, or reverse charge rechargeable batteries.
- Never open, dismantle or damage rechargeable batteries.
- Never throw rechargeable batteries in fire.
- Never expose rechargeable batteries to excessive heat or direct sunlight.
- If the device becomes unusually hot, emits an odour, changes colour, releases gas or smoke, or if the charging process takes longer than normal, stop using and charging the device, and contact our customer service team.
- If the rechargeable battery is damaged or leaking, avoid contact with skin, eyes, or mucous membranes. In the event this occurs, immediately flush the affected areas with plenty of clean water and seek medical help.

Danger for children, persons with impaired physical, sensory, or mental capacities

e.g. partially disabled persons, older persons with reduced physical and mental capacities, and those lacking experience and knowledge (e.g. older children).

- The device may be used by children aged 8 years and above, persons with impaired physical, sensory, or mental capacities, and those lacking experience and knowledge if they are supervised or have been instructed in how to safely use the device and have understood the risks associated.
- Cleaning and user maintenance must not be performed by children.
- Keep children under 8 years of age away from the device and all its accessories.
- Do not leave children unattended while using the device.
- Do not allow children to play with the device or any of its accessories.

Risk of injury

Handling or using the device improperly may result in injury.

- Never direct the water jet under the tongue, into the throat, eyes, ears, nose, or other delicate areas. The device produces pressures that may cause damage in these areas.
- Remove any oral jewellery prior to using the device.
- When you first use the device, you may experience moderate discomfort or bleeding, this is normal and will usually pass after one weeks of daily use. Try using a lower pressure setting until you become more accustomed to the device.
- Do not press the floss tip forcefully on the gums or teeth.
- If you have any oral conditions, please consult your dentist before using the device.
- If you are due to undergo oral surgery or have undergone any oral surgery in the last three months, consult your dentist before using the device.
- Do not use the device if you have an open wound on your tongue or in your mouth.
- Stop using the device and seek advice from your dentist or doctor if you experience: excessive bleeding after use, bleeding that persists beyond one week of use, or excessive or prolonged discomfort or pain.
- Stop using and replace the floss tip if it becomes worn or damaged, to prevent injury to the mouth and tooth area. We recommend replacing floss tips every 6 months or sooner if signs of wear appear.
- Do not operate the device without a floss tip installed.
- Lay the USB cable so that it does not pose a tripping hazard.
- Do not operate the device where aerosol (spray) products are being used or where oxygen

is being administered.

Risk of damage

Improper handling of the device may damage it.

- Place the device on an easily accessible, level, and dry surface.
- Keep the device and all its accessories away from hot surfaces and open flames.
- Be careful not to drop the device. If this does happen, do not use it anymore; instead, contact our customer service team. The floss tip must be replaced, even if no damage is visible.
- The device shall only be used with cold or warm water (below 40 °C), diluted Ordo mouthwash, or other dental professional recommended solutions. Do not use pure mouthwash with the device, only diluted mouthwash is permitted. After using diluted mouthwash, run a full tank of clean water through the device to clean it.
- Do not use iodine, salt solution or water insoluble concentrated essential oils in the device. Use of these can reduce device performance and will shorten the life of the device.
- Do not use the device while it is immersed in water.
- Charge the device in an upright position and with the water tank attached.
- While the water tank is safe to clean in a dishwasher, do not clean any other part of the device in a dishwasher.
- To prevent scratching the main body of the device, use a soft, damp cloth for cleaning.
- The device should only be cleaned with clean water. Do not use abrasive or solvent-based cleaning products.
- If you will not use the device for a long time, empty the water tank, clean and dry all parts, and store in a cool and dry place.
- Do not use the USB cable with any other device, as this may cause damage.

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