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

ELECTRONIC PERSONAL SCALE VT-1951

This model of scale uses BIA bioelectric body total resistance (impedance) analysis, which provides highly accurate calculations of the percentage of fat tissue, water, bone and muscle tissue to calculate individual body mass, necessary daily calorie intake, based on user's personal data (sex, height, age) input into the unit memory.

DESCRIPTION

1. Weighing platform
2. Contact areas
3. Weight measurement unit selection button «Unit»
4. Data selection buttons «▲▼»
5. Data selection and confirmation button «SET»
6. Display
7. Battery compartment lid

Display symbols

8. Digital indications
9. User's height in centimeters «cm»
10. Symbol of the sex: female «♀»
11. Symbol of the sex: male «♂»
12. Weight symbol in kilograms «kg»
13. Percentage symbol «%»
14. Weight symbol in stones and pounds «st:lb»
15. Body mass index symbol «BMI»
16. Symbol of measurement of the percentage of adipose tissue «obese» «**OBESE**»
17. Symbol of measurement of the percentage of adipose tissue «above normal» «**OVERFAT**»
18. Symbol of measurement of the percentage of adipose tissue «norm» «**HEALTHY**»
19. Symbol of measurement of the percentage of adipose tissue «below normal» «**UNDER FAT**»
20. Fat percentage measurement symbol «**FAT**»
21. Water percentage symbol «**TBW**»
22. Bone mass symbol «**B**»
23. Muscle bulk symbol «**M**»
24. Recommended daily calorie intake «**KCAL**»
25. User's age symbol «**Age**»

SAFETY PRECAUTIONS

ATTENTION!

- People having pacemakes and other medical devices should not use the scale.
- Pregnant women should not use the scale.

SAFETY MEASURES AND OPERATION RECOMMENDATIONS

Before using the personal electronic scale (hereinafter - the scale), read this instruction carefully and keep it for further reference.

Use the scale for its intended purpose only, as specified in this manual. Mishandling the scale can lead to its breakage and cause harm to the user or damage to their property.

- Handle the scale with care, as any measuring instrument.
- Avoid exposing the scale to direct sunlight and do not drop the scale.
- Do not expose the scale to high or low temperatures.
- Do not keep or use the scale in a room with high humidity (above 80%), make sure that the scale body does not come into contact with water or any other liquids as this may cause scale reading malfunction or its damage.
- Use the scale away from heating units.
- Place the scale on a flat dry and steady surface.
- For more precise readings do not use the scale on carpets or rugs.
- Do not step on the scale with wet feet or if the weighing platform is damp. You can slip and get injured.
- Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and only then step on it with the both feet.

- Stand still while weighing. Do not jump on the weighing platform.
- If you are not planning to use the scale for a long time, remove the batteries from the battery compartment.
- Strictly follow the polarity according to the markings when replacing the batteries.
- Replace the batteries in time.
- For children safety reasons do not leave polyethylene bags, used as a packaging, unattended.

Attention! Do not allow children to play with polyethylene bags or packaging film. **Danger of suffocation!**

- Do not leave children unattended to prevent using the unit as a toy.
- This unit is not intended for usage by children. Place the unit out of reach of children.
- The unit is not intended for usage by physically, sensory or mentally disabled persons (including children) or by persons lacking experience or knowledge if they are not under supervision of a person who is responsible for their safety or if they are not instructed by this person on the usage of the unit.
- Do not exceed the maximum capacity of the scale which is 180 kg.
- Do not repair the unit by yourself. Do not disassemble the unit by yourself, if any malfunctions are detected, remove the batteries from the battery compartment and apply to any authorized service center at the contact addresses given in the warranty certificate and on the website www.vitek.ru.
- Transport the unit in the original packaging only.
- Keep the unit out of reach of children and disabled persons.

THIS UNIT IS INTENDED FOR HOUSEHOLD USE ONLY. ITS COMMERCIAL USAGE AND USAGE IN PRODUCTION AREAS AND WORK SPACES IS PROHIBITED.

IMPORTANT

Electromagnetic compatibility

Electronic scale can be sensitive to the electromagnetic emission of other units located in close proximity (such as cell phones, portable radio transmitters, remote controllers and microwave ovens). Should the signs of such emission appear (false or inconsistent data readings on the display), use the unit away from the source of interferences.

BEFORE THE FIRST USE

After unit transportation or storage at low temperature keep it for at least 3 hours at room temperature before switching on.

- Wipe the scale with a soft, slightly damp cloth and then wipe it dry.

BATTERIES INSTALLATION AND REPLACEMENT

ATTENTION!

Batteries leakage can cause injuries or damage to the unit. Please follow the instructions to avoid damage:

- use the «AAA» type batteries only;
- do not use old and new batteries at the same time or different types of batteries;
- do not recharge the batteries;
- observe the polarity strictly while installing the batteries;
- remove the batteries if you are not planning to use the unit for a long time;
- apply to specialized centers for further recycling of the batteries;
- avoid short circuiting the battery clips.
- Remove the battery compartment lid (7) and insert 2 «AAA» batteries (supplied with the unit), strictly following the polarity.
- Install the battery compartment lid (7) back to its place.

ENGLISH

Battery replacement

- When the batteries are low, the symbol «**LO**» will appear on the display (6).
- Open the battery compartment lid (7), replace the «AAA» batteries with new ones, strictly following the polarity, and close the battery compartment lid (7).
- If you are not planning to use the scale for a long time, remove the batteries from the battery compartment.

RECOMMENDATIONS

- Try to install and use the scale in the same place.
- Weigh yourself at the one and the same time of the day.
- We recommend measuring the percentage of fat and water in the evening before evening meal.
- 2-3 hours should pass after last meal.
- At least 10-12 hours should pass after intensive physical exercises.
- Percentage of fat in the body depends on different factors and can change during the day.
- When calculating percentage of fat and water of children younger than 10 years the results may be inaccurate.
- Professional sportsmen, people inclined to edemas and people with nonstandard body proportions can show inaccurate results of calculation of fat and water percentage.

Note: Fat and water percentage data can be calculated only if personal data (sex, age, height) are entered correctly and if you step on the scale barefeet.

WEIGHT MEASUREMENT

Weight measurement unit selection

- To select the weight measurement unit, consequently press the button (3) «**Unit**» on the bottom part of the scale body.
- Measurement units are shown on the display (6) with the symbols:
«**kg**» – weight measurement in kilograms (12);
«**lb**» – weight in pounds (14);
«**st**» – weight in stones (14).
- Place the scale on a flat solid surface; if the floor is uneven or covered with carpet or rug, the weighing accuracy may decrease.
- Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and only then step on it with the both feet.
- Step on the scale and dispose your weight evenly. Stand still while being weighted. Do not jump on the weighing platform. The digital indications (8) on the display (6) will show your weight in the selected measurement units.

ATTENTION!

- If a symbol «**Err**» appeared on the display (6), it means that the scale is overloaded.
- Immediately step off the weighing platform (1) to avoid damaging the scale.
- The scale will be automatically switched off in 42 seconds after switching on.

ENTERING PERSONAL DATA

- Switch the scale on by pressing on the weighing platform (1). Press the «**SET**» button (5), the flashing symbols (8) on the display (6) will show one of the user numbers «**P1**», «**P2**», ..., «**P13**».
- Select the desired user number by pressing the buttons (4) «**▲▼**». You can enter and save data for 13 users. Confirm your selection of the user number by pressing the «**SET**» button (5).
- The flashing symbol of the user's sex will be shown on display (6). By pressing the buttons (4) «**▲▼**» select the symbol (10) «**♀**» for woman or the symbol (11) «**♂**» for man, press the «**SET**» button (5) to confirm.

- The flashing digital indications of user's height will be shown on the display (6), the digital values of height are entered in values corresponding to the selected weight measurement units «**kg/st/lb**». If the weight measurement units (12) are set in «**kg**», the user's height values are entered in centimeters (9) «**cm**». Set the height values by pressing the buttons (4) «**▲▼**», press the «**SET**» button (5) to confirm.
- The flashing values of the user's age (25) «**Age**» will appear on the display (6). Set the user's age by pressing the buttons (4) «**▲▼**», press the «**SET**» button (5) to confirm.

WEIGHT MEASUREMENT WITH CALCULATION

RESULTS

- Switch the scale on by pressing on the weighing platform (1). Press the «**SET**» button (5), the flashing symbols (8) on the display (6) will show one of the user numbers «**P1**», «**P2**», ..., «**P13**».
- Select the desired user number by pressing the buttons (4) «**▲▼**». The user's number symbol will appear on the display (6).
- Wait until the symbols «**0.0 kg**», «**0.0 lb**» or «**0:0 st**» appear on the display (6), depending on the selected measurement units.
- Step with your bare feet on the weighing platform (1) so that your feet are on the contact areas (2), stay still during the weight measurement and wait until the weight indications appear on the display (6).

Note: If you don't step on the scale within 42 seconds, it will switch off automatically.

- Keep standing still on the weighing platform (1), the display (6) will show moving symbols «**0...0**», while the data is being calculated.
- The calculated data will be shown on the display (6) in the following order:
 - «**FAT**» - fat tissue percentage (20);
 - «**TBW**» - water percentage (21);
 - «**MB**» - muscle bulk percentage (23);
 - «**BM**» - bone mass percentage (22);
 - «**KCAL**» - recommended daily calorie intake (24);
 - «**BMI**» - body mass index (15);
 - after the calculated results, your weight will be shown once again;
 - after that the scale will switch off.

As a result of calculations, the display (6) can show the following symbols:

- obesity symbol «**OBESE**» (16);
- overnutrition symbol «**OVERFAT**» (17);
- healthy body symbol «**HEALTHY**» (18);
- excessive fat symbol «**UNDER FAT**» (19);
- the appearance of the symbol «**ErrL**» indicates an error in the percentage of body fat, or very low body fat.
- the appearance of the symbol «**ErrH**» indicates an error in the percentage of body fat, or very high body fat.

CALCULATION RESULTS

- «**FAT**» - the result of calculating the percentage of fat tissue in the body, this parameter is not constant and may change during the day. The percentage of fat tissue in the body depends on lifestyle, physical activity, general health and the amount of fluid in the body.
- When analyzing calculated readings of fat tissue in the body, consider the following factors:
 - the fat percentage indications may be higher after having a meal or a drink,
 - the fat percentage indications may be less after taking a bath or a shower,
 - if the fluid content in the body decreases, then the fat tissue percentage indications will change.

Fat tissue correlation table

Age	Women (%)				Men (%)			
	Under Fat	Healthy	Overfat	Obese	Under Fat	Healthy	Overfat	Obese
< 20	< 23	23-33	33-42.3	> 42.3	< 18	18-28.2	28.3-35.6	> 35.6
20-30	< 21.6	21.6-33	33-41.2	> 41.2	< 16.9	16.9-27.2	27.2-34.6	> 34.6
30-40	< 20.5	20.5-31.9	31.9-40	> 40	< 15.9	15.9-26.2	26.2-33.7	> 33.7
40-50	< 19.4	19.4-30.8	30.8-39	> 39	< 14.9	14.9-26.2	26.2-32.7	> 32.7
50-60	< 18.3	18.3-29.7	29.7-37.8	> 37.8	< 13.9	13.9-24.3	24.3-31.7	> 31.7
> 60	< 17.2	17.2-28.6	28.6-36.7	> 36.7	< 12.9	12.9-23.3	23.3-30.7	> 30.7

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

- «**TBW**» - water percentage in a human body is the percentage of liquid in the human body to total weight. The human body consists of 50-65% water, so water plays a primary role for tissues and organs.

This scale helps to control the water percentage in a body.

Water correlation table (%)

Age	Women (%)			Men (%)		
	Low	Normal	High	Low	Normal	High
10-15	< 57	57-67	> 67	< 58	58-72	> 72
16-30	< 47	47-57	> 67	< 53	53-67	> 67
31-60	< 42	42-52	> 67	< 47	47-61	> 61
61-80	< 37	37-47	> 67	< 42	42-56	> 56

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

- «**BM**» - muscle tissue percentage depends on the sex, weight, height, age of the person and on their lifestyle.

Average values for women are: - > 34%, for men: - > 40%

- «**B**» - bone tissue mass.

The human skeleton consists of hard bone tissues that protect the internal organs, and together with muscle tissue provide a person's ability to move.

The calculation of bone mass is based on the principle of determining the mass of mineral substances, in particular calcium and other mineral substances. Research results have shown that development of muscle tissues makes the human skeleton stronger. Take this into account when calculating physical activity and diet. The mass of bone tissue practically does not change, but its slight decrease is possible under the influence of many factors (sex, weight, height, age).

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Recommended daily calorie intake table (%)

Age/Weight		< 17 years	18-29 years	30-50 years	> 50 years
Men	< 50 kg	1350	1200	1115	1075
	50-65 kg	1485	1320	1227	1183
	65-80 kg	1890	1680	1561	1505
	80-95 kg	2295	2040	1896	1828
	> 95 kg	2700	2400	2230	2150
Women	< 50 kg	1265	1180	1085	1035
	50-65 kg	1392	1298	1194	1139
	65-80 kg	1771	1652	1519	1449
	80-95 kg	2151	2006	1845	1760
	> 50 kg	2530	2360	2170	2070

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

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The body mass index allows to assess the degree of conformity of a person's weight with their height and thus indirectly tell whether the person's body weight is insufficient, normal or excessive.

Apply the body weight index with caution, only for a rough estimate – for example, an attempt to assess by its means the body build of professional athletes can lead to incorrect results (high index value in this case is accounted for by the athlete's well developed musculature).

In the table below the body mass index data are given according to the WHO estimation

«BM1» parameter	Under Fat	Healthy	Overfat	Obese
Body weight index	< 18.5	18.5-25	25-30	> 30

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

COMPARISON OF THE RESULTS OF WEIGHING AND CALCULATION OF PARAMETERS

- The results of parameter calculations and simple weighing are only relevant for long-term comparison of results, since short-term changes in body mass index are usually caused by a slight loss of fluid in your body tissues.
- Evaluate results based on changes in total body weight, percentage of body fat, and fluid content in the tissues, during the period of time when these changes occurred. Short-term indication changes during a day can be significantly different from the indication changes within a week or a month.
- Short-term weight loss and simultaneous growth or constant fat percentage in body tissues may be caused by decreasing of liquid amount as a result of physical exercises, visiting sauna or quick weight loss.
- Medium-term increase of weight, decrease of fat percentage or constant fat tissue amount in the human body may be caused by increasing of muscle bulk.
- Simultaneous decrease of total weight and fat tissue amount shows that your diet is working.
- Combine diet with different physical exercises. In this case you will observe the medium-term increase of muscle mass in body tissues.

CLEANING AND CARE

- Clean the weighing platform (1) with a soft slightly damp cloth, and then wipe dry.
- Do not use abrasive cleaners or solvents to clean the scale.
- Do not immerse the scale into water or any other liquids.
- Keep the scale in a dry cool place out of reach of children and disabled persons.

ENGLISH

DELIVERY SET

Scale – 1 pc.
«AAA» battery – 2 pcs.
Instruction manual – 1 pc.
Warranty certificate – 1 pc.

TECHNICAL SPECIFICATIONS

Power supply: 3 V, 2 «AAA» batteries
(included in the delivery set)
Maximum weight limit: 180 kg
Minimum weight limit: 5 kg.
Division scale*: 100 g

* For household scales the extent of error is not specified, so it is not stated in technical specifications.
The term «margin of error» is used for industrial, analytical, and medical scales. These scales have pattern approval certificates which confirm the measurement accuracy. Moreover, specially appointed organizations perform annual "measurement accuracy confirmation" (checking) to confirm the scales accuracy.

RECYCLING

For environment protection do not throw out the unit and the batteries (if included), do not discard the unit and the batter-

ies with usual household waste after the service life expiration; apply to specialized centers for further recycling.

The waste generated during the disposal of the unit is subject to mandatory collection and consequent disposal in the prescribed manner.

For further information about recycling of this product apply to a local municipal administration, a disposal service or to the shop where you purchased this product.

The manufacturer preserves the right to change design, structure and specifications not affecting general principles of the unit operation without a preliminary notification due to which insignificant differences between the manual and product may be observed. If the user reveals such differences, please report them via e-mail info@vitek.ru for receipt of an updated manual.

Unit operating life is 3 years

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC Directive 2014/30/EU and to the Low Voltage Directive 2014/35/EU.

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