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MANUAL
ELECTRIC WATER FLOSSER
SD-WF22503PP

ELECTRIC WATER FLOSSER

Thank you for purchasing a product manufactured with Swiss Diamond standards. Our products are made in accordance with the highest demands of a quality, utility and design. We hope that you will continue to choose our products in the future.

Before using the product, make sure to read this manual carefully, it contains an important safety information and recommendations for use and correct maintenance.

Please save the instruction manual with a warranty card, cash receipt and if possible, with a gift box.

Introduction:

1. Description
2. Water flosser head types
3. Important safeguards
4. Instructions for use
5. Features and properties
6. Charging of the battery
7. Cleaning and maintenance
8. Specification
9. Warranty
10. Certification

1. DESCRIPTION

Standard head

Place for heads

On/Off button

Mode selection button

Panel displaying selected mode

Battery power indicator

Charging port

Removable water tank

Release button

Suction tube

Water tank lid

Anti-slip feet

2. WATER FLOSSER HEAD TYPES

CLASSIC HEAD

It is used for daily cleaning of teeth and interdental spaces from bacteria and food remains. Also suitable for removing food remains from imitative constructions.

2 pcs in set.

ORTHODONTIC HEAD

It is suitable for places hardest to reach: braces, crowns, dental bridges and implants. During use hold the head between your teeth for a while to clean the area between them and around the bracket (crown). It's recommended to use this head with a small water pressure so as not to damage the braces or other systems.

1 pc in set.

CAUTION! The water coming from the head should not be too hot or cold. It is recommended to use water at room temperature.

PERIODONTAL HEAD

It is used for cleaning of gum line and gum pockets. Place the soft silicone tip of the head against the tooth at 45C angle and gently slide the tip under the gum line into the gum pocket.

1 pc in set.

TONGUE HEAD

It is used for daily cleaning of tongue. It removes the bacterial plaque and freshens breath. Place the head on the tongue, turn on the flosser and move the head to the tip of tongue with a little effort.

1 pc in set.

CAUTION! It is recommended to replace the heads every 6 months with regular use.

3. IMPORTANT SAFEGUARDS

Read all instructions and important safeguards to avoid damage or injury when using this electric water flosser:

1. The device should not be used in any way other than the one specified in the manual. The water flosser is made only for oral care: cleaning areas between the teeth, gums and imitative constructions (braces, crowns, etc.).
2. Before plugging the device, make sure that the voltage shown on the device is the same as the one in the outlet.
3. If the device is damaged in any way or makes unusual noises after being powered on (in the area of head, main body or charger), do not use it. It is forbidden to try to repair it on your own. The device does not contain any elements made for the user's repairs.
4. The device should not be used by the people (including children) with physical, sensory or mental impairments and those with absence of knowledge or experience necessary to operate such devices, unless they are being supervised or instructed about usage of device by the person responsible for their safety.
5. The children should be instructed to not use the device as a toy.
6. While unplugging the device, always grab and pull the plug, not the cable.
7. The charger and power cord should be used and stored away from the heat sources, water and damp.
8. The manufacturer does not take responsibility for the damages caused by an improper use of the product.
9. Do not use the heads other than the ones recommended by the manufacturer.
10. Do not use the charger while in bath or under the shower.
11. Do not use this device in places where it can land or be

thrown into the water. If the water flosser has been dropped on the floor, replace the head before next use, even if there is no visible damage.

12. Do not submerge the handle or the charger in water or any other liquid. After cleaning and before plugging the device, make sure that every component is dry.

13. The device is for the home use only. It is not made for professional use.

14. The device is for oral care only. Do not use for other purposes. Do not chew or gnaw on the heads.

15. If you are undergoing treatment for any oral cavity disease or had jaw/gum surgery within the last 2 months, make sure to contact your dentist before using the device.

16. Consult your dentist if your gums start bleeding or you feel pain after using the water flosser.

17. If the user has a pacemaker or similar device, consult a dentist or product manufacturer before using the device.

18. Replace the heads every 6 months or sooner (if you notice the signs of wearing out).

19. Do not clean the heads in the dishwasher.

CAUTION! The device has a waterproof handle and is protected from moisture, but it is not recommended to use it under water!

20. Do not direct water under the tongue, in the ear, nose, eye or other sensitive areas. This device delivers water under pressure which can cause serious damage to these areas. Also, it is not recommended to use the device if you have wounds, stomatitis on the tongue.

21. Do not switch the device if the water tank is empty.

22. Fill the water tank by clean drinking water or solutions recommended by your dentist and nothing else. The liquid mustn't be with sediment, viscous, otherwise the device will become clogged.

23. Do not use solutions that contain iodine and bleaching elements.

24. Before using the device, remove all foreign objects from the oral cavity (piercing, jewelry, etc.).

25. Pour out any remaining water after using the device and dry the water tank before joining to the main body.

26. Water temperature should be not higher than 40C° to avoid burns of skin and mouth.

4. INSTRUCTIONS FOR USE


01 Assembly

Before using the device, install the necessary head on top of the handle until it clicks.


02 Removal, installation and filling the water tank

Holding the body of the device with one hand and the water tank with other hand, slightly turn the body of the device counterclockwise and disconnect the water tank (see picture below and on the device).

Open the water tank lid and pour water in water tank as shown on the picture. After filling the water tank, close the lid tightly and connect water tank to the main body. You may also fill the water tank without disconnecting it from the main body.


03 Switching on the device

Press the On/Off button to switch the device on.

04 Mode selection

To choose the operating mode, press the mode selection button until the corresponding indicator light on the panel lights up.

Water flosser offers 5 working modes, allowing you to adjust it to your needs:

Strong: for deep cleaning of the oral cavity from food remains and plaque. Suitable for people with healthy gums. In this mode water is delivered with the maximum pressure.

Normal: for daily cleaning. Suitable for most people. In this mode water is delivered with the medium pressure.

Soft: for gentle cleaning. Suitable for sensitive teeth and gums.

Pulse: for cleaning and massaging gum pockets. Gum massage enhances blood microcirculation which is a good prevention of periodontitis and periodontal disease.

Short: for the prevention of calculus plaque, gum massage and point impact on the problem area. In this mode water is delivered at short intervals.

Advice: it is recommended to choose the soft mode during your first usage.

05 Cleaning

Bend over the sink and hold the device strictly in a vertical position, aiming the head at your teeth and gums. Do not press On/Off button until the head is in your mouth. Start cleaning from the back of the jaw and slowly move along the teeth line, stopping at the places between teeth for better disinfection. Cover your mouth slightly during irrigation to prevent splashing.

5. FEATURES AND PROPERTIES

Waterproof level: IPX7

High waterproof level for comfortable and safe use.

Remembering the last mode chosen

After turning the device on again, it will be set to the previously chosen cleaning mode.

Two-minute timer

Water flosser can turn itself off automatically after 2 mins of usage.


Up to 50 minutes of work

After a single charge, the water flosser can be used for 14 days (work in Strong mode).

6. CHARGING OF THE BATTERY

1. Before using the water flosser for the first time, make sure that the device is **fully** charged.
2. Plug in the charger to the charging port of the water flosser.

3. During the charging procedure, the battery indicator should **blink** in red color.
4. The indicator will **continue blinking** until the water flosser is **fully** charged.
5. Once charged, the device will indicate the **full** battery by **constant white** light.
6. A **full** charging cycle lasts from 4 to 5 hours, **allowing** for 14 days of continuous work (depending on the selected mode).
7. Rapid **blinking** of the battery indicator means that the battery power is running **low**.

CAUTION! The device can only be used when it is charged! Mains usage (during charging) is impossible.

7. CLEANING AND MAINTENANCE

1. Clean the water flosser regularly. Wash the main body, water tank and heads with running water after each use.
2. Pour out any remaining water and carefully dry the water tank before assembling.
3. Leave the water tank lid open.
4. Store your device in a cool and dry place

Environment friendly disposal

Please remember to respect the local regulations: hand in the non-working electrical equipment to an appropriate waste disposal center.

8. SPECIFICATION

Product type:	water flosser SWISS DIAMOND
Working modes:	strong, normal, soft, pulse, short
Switch type:	a button
Run time:	50 min
Charging time:	4-5 hours
Power:	8 W
Battery:	Li-ion, 1500 mA h
Battery indicator:	flashing light
Noise level:	≤75dB
Waterproof level:	IPX7

9. WARRANTY

- 36 months of warranty is valid for this device.
- During the warranty period, the manufacturer is obliged to eliminate any manufacturing defects by repairing, replacing parts or replacing the entire product.
- The warranty comes into force only if the date of purchase is confirmed by the store's seal on the warranty card.
- This warranty is valid if the product has been used in accordance with the instruction manual, has not been repaired, disassembled or damaged as a result of improper usage.
- This warranty does not cover to body trim details, natural wear, consumables and accessories (heads, protective caps, etc.).

The full list of service centers you may find on the web-site: www.inhousewares.com

Life time is 5 years.

The month and year of manufacture are indicated on the bottom of the device and packaging.

The manufacturer has the right to make changes in the design, configuration, packaging, and also in the technical features of the product in order to improve products without further notice of these changes.

10. CERTIFICATION

This product meet the requirements of the technical regulations of the Customs Union, all required European and Russian safety and EMC standards.

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