Smartwave[®]

Air compression
therapy device
Model: 600

Smartwave[®]

Address: B1/1-5F, Tonglu Industrial Area, 70# Tongxin Road, Tongxin Community, Longgang District, Shenzhen 518116, Guangdong, China

www.smartwave-europe.eu

User manual

Thanks for purchasing our product.
Please read the instructions carefully prior to the use
and keep it well for further reference.

EN 2-15

RU 16-29

Intended User Population

This device is intended for use by adults.

Intended User

The PATIENT is an intended OPERATOR Applied Part: Sleeves.

Indication for Use/intended use Smartwave 600

- relaxing massage
- removal of edema of all kinds
- treatment and prevention of varicose veins
- stimulation of blood circulation
- stimulation of lymph flow
- venous insufficiency
- body shaping
- prevention and elimination of cellulite
- getting rid of excess weight
- increased skin tone
- recovery after exercise
- recovery from injuries

Intended User

Do not use this product if you are experiencing inflammation, an infection, pain of unknown origin, bleeding (internal or external) at or near the site of application, or if you have a wound at or near the site of application. Do not use this product if you are under the care of a physician or have a condition requiring the use of any medical device.

- Do not use this product on sensitive skin.
- Do not use this product if you have any of the following conditions:
- Acute pulmonary edema
- Acute thrombophlebitis
- Acute congestive cardiac failure
- Acute infections
- Deep vein thrombosis (DVT)
- Wounds, lesions, or tumors at or near the site of application
- Bone fractures or dislocations at or near the site of application
- Presence of the arterial insufficiency

Side effect

You may experience discomfort when the air pressure is too high, in this case decrease the pressure, procedure must be comfortable.

Precautions


It is important that you read all the warning and precautions included in this manual because they are intended to keep you safe, prevent injury.

Consult with doctor if a user is as below:

- A person who is in pulmonary edema, heart disease, vein cancer, high blood pressure, high fever.
- A person who has skin disease such as skin burn, dermatitis, bedsore, purulent disease, malignant tumor and etc.
- A person who took skin surgery (skin graft) on or near the treatment area.
- A person who has blood vessel disease such as series artery hardening or other angina, cardiac infarction and etc.
- In the case of the suspected presence of blood clot.
- The treatment area with physical deformity or pins inserted, (artificial joint, beads, metal, and etc.)
- Just experienced the varicose veins operation.
- Recovering patients since after the surgery, pregnant woman, children.
- In the case of pain increasing after use of this machine.

WARNINGS

	Please avoid the same time with other medical equipment and medical treatment implants, such as pacemaker; artificial heart and etc., for safety. Users with too high blood pressure or obstacle in blood circulation are not allowed to use.
	Do not place heavy things on the power cord. Also, do not press it with a table or a chair. It may cause an electric shock, an electric leakage, or a fire.
	When plugging or unplugging the power cord from the wall outlet, please remove it with dry hands. Otherwise it may cause electric shock.
	Please wear something over when using these sleeves, to avoid a direct contact between the skin and the sleeves.
	Be careful not to drop the device on the floor because much vibration may cause functional damage to it.

	Do not use it in an environment with too much electromagnetic wave. It may cause malfunctioning to the device and damage to the user.
	Be careful that water or other substances cannot be entered into the device, do not immerse the device in water.
	While using this device, do not attempt to walk or move around.
	If you have not used the device for a long time, or if you wish to use it in cold areas. The device may not operate normally. In such case, use the device after operating it three or four times with different pressure.
	Please strictly follow the cleaning instructions as explained in this manual to clean the device, otherwise it may cause the product external damage such as discoloration.
	Please plug the device to a exclusive power outlet, do not plug other devices to a same outlet for use at once due to the risk of a fire, electric shock or other injuries.
	If this device malfunctions or fails to operate, do not attempt to open, repair, or disassemble the main unit or this device due to the risk of a fire, electric shock, or injuries to the user.
	Using this device on the lowest settings in first use, and increase the treatment time and pressure slowly according to your needs and feeling...
	Do not get the air hose twisted or bent.
	Do not use this device if you are experiencing inflammation, an infection, pain of unknown origin, or bleeding (internal or external) at near the treatment position.
	Do not machine wash the sleeves or submerge them in water. Follow instructions in this manual to clean them.
	Stop using the device immediately when you feel pain or uncomfortable. Be careful when you zip up the sleeve, as clothing, hair or skin may get stuck in the zipper.
	Check up the suitable voltage for the use of this device. Otherwise it may cause a function failure, electric shock, fire, or injuries

	Do not use by purpose except for massage. Keep it away from infants and children due to the risk of electric shock or other injuries.
	This device has not been tested for use in an MR environment and should not be used exposed to MR environments while patients are wearing the device. Keep it outside the MRI scanner room. Use the device according user manual's instruction.

⚠ WARNINGS

- The equipment are not serviced or maintained while in use with the patient.
- Prohibition of Sitting and Stepping on the device, Prohibition of Pushing or Dis-assembly of the device.
- For replacement accessories, please note that only the original parts from the manufacturer should be used.
- Mains plug is used as the isolation means from mains supply. Do not position the device to make it difficult to disconnect mains plug from power socket.
- Avoid strangulation resulting from power supply cord and joint air hose.
- Keep the device away from pets, children to avoid unintentional destruction.
- Keep the device well storage to avoid accelerated ageing caused by lint, dust and sunlight -Failure of temperature sensor will lead to the shell over heat.
- Failure of pressure sensor will lead to the air pressure increases .

WARNING: Using microwave oven, HF surgical equipment, magnetic resonance imaging or other radio radiant equipment near this product may cause malfunction or lead to loss of essential performance.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the air compression therapy device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

List of cables and maximum length of cables is as follows:

Cables name	Cable length	Whether shielding
Power cord	1.8m	No

Caution:

Security, anti-theft, and radiofrequency identification (RFID) devices. Some electromagnetic anti-theft systems and metal detectors such as those used at entrances or exits of department stores, libraries, and other public places, and airport security screening devices may affect the Air compression device. Additionally, RFID devices, which are often used to read identification badges, as well as some tag deactivation devices, such as those used at payment counters at stores and loan desks at libraries, may also affect the Air compression therapy device. Please do not use Air compression therapy device near these places.


If you experience severe pain or any unusual symptoms:

- Stop the device by pressing the standby button.
- Disconnect the connector of the air Hose from either the main Unit or leg cuffs.
- Remove the cuffs from your limb(s).
- Promptly consult your licensed healthcare practitioner, as required.

Product description


The product is composed of main unit, leg cuffs and joint air hoses.

Main unit


No.	Name	Description
1	Control panel	User interface, for the user to operate the equipment; display treatment information.
2	Air outlet	Used to connect the Hose connector.
3	Power plug	Used to connect an AC socket.
4	Power cord	A wire that transmits current to main unit
5	Carry handle	Handle, convenient for users to move the device.

Accessories

No.	Name	Instruction
1	Attachment connector	Used to connect to the Junction box of Joint air hose
2	Zipper	Zip up before use so that cuff can wrap the legs
3	Air Chambers	The part where the air is inflated and deflated, 6 air chambers for each leg cuff.
4	Junction box	Used to connect to the Attachment connector of Leg cuff.
5	Air hose	The hose that transmits air to air chambers.
6	Hose connector	Used to connect to the Air outlet of main unit.

Control panel of main unit

No.	Name	Instruction
1	LED screen	• Display setting mode: When press "Mode" button, It will display the setting Mode (P1 or P2 or P3) for 2s; • Display pressure value: 2s after pressing the mode key, it will automatically switch to display the current air pressure value(50~199 mmHg).
2	Time button	Choose the massage timing: 10/20/30 minutes can be selected. Each timing has its own indicator (10Min/20Min/30Min indicator).
3	Standby Button	Turn on/off the main unit
4	Chamber 1~Chamber 6 button	Treatment site control button (Chamber 1, Chamber 2, Chamber 3, Chamber 4, Chamber 5, Chamber 6 button control the corresponding treatment site from foot to thigh)
5	Mode button	To select a mode.
6	Pressure knob	Adjust the pressure level: rotate to the left, pressure decreases, rotate to the right, pressure increases
7	Diagram of treatment sites	Six blue indicators, indicating the corresponding treatment position.


A General Specifications of device

Product name	Air Compression Therapy Device	Model	600
Dimension (main unit)	258 x 192 x 119mm	Weight (with all accessories)	3.5KG
Pressure Range	50~199 mmHg	Power supply	AC Input: 100~240 Volts 50/60HZ, 1.0A Max.
Power Consumption	35w	Number of chambers	6 chambers
Mode of operation	Continuous Operation	Mode	3 modes
Against electric shock	Applied part	Treatment Position	Feet, calves, knee and thigh
Treatment time	10/20/30 minutes	Shell materials	ABS plastic
Cuff size	1015 x 380mm		
Operating Conditions	Temperature: +5° C to 40° C; Humidity : 5% to 90% non-condensing		
Storage Conditions	Temperature: -20° C to 55° C; Humidity : 5% to 90% non-condensing. Keep dry and avoid direct sunlight exposure.		

How to use


Setup the system

- step 1** Pull down the zipper of cuff, then put your whole leg in and wrap your legs well by zipping it up. There is no left-right difference between the sleeves.

step 2 Plug "attachment connector" of Leg Cuffs to the "junction box connector" of Joint Air Hose completely; then properly plug the "hose connector" of the joint air hose to the "Air outlet" of main unit.

Note: You should hear a light click when you completely plug the hose connector into the main unit

step 3 Connect the power plug of the main unit into an AC socket (the mains supply).

step 4 Press the "Standby" button of the main unit to start massage.


It will begin with Sequence mode and the air pressure of 50mmhg (The lowest), 10-minute timing by default. It will shut off automatically after 10 minutes treatment if no timing selected (can manually re-start the device).

Select the mode


Press "Mode" button to switch and enjoy the different massage modes, the LED screen will show the select mode (Sequence mode: display P1, Circulation mode: display P2, Combination mode: display P3) for 2s after press the mode button, and the corresponding LED indicator will flash. Adjusting mode during operation will not make the timing re-start.

The massage sequence of each mode is as follows:


Mode1: Sequence Mode

Mode 2: Circulation Mode

Mode 3: Combination Mode

Select the treatment area

Press "Chamber 1", "Chamber 2", "Chamber 3", "Chamber 4", "Chamber 5" or "Chamber 6" button to stop or restart the massage of its matched chambers.

The indicator will flash when its matched chamber is working. When the chamber stops working, its indicator stays off;

All chambers work by default when the device is turned on, all indicators of treatment areas show a fixed blue light when the main unit is turned on; Adjusting the treatment area during treatment does not result in timing re-start. Adjusting the treatment area will not disturb the sequence of mode, the treatment area that are not selected will be skipped in sequence working process.

Note: The main unit will shut off automatically after 1 minute when all chambers are turned off.

Adjust the intensity and time

Time adjustment:

Press "Time" button to choose a operating timing from 10/20/30 minutes, the LED indicator of the timing applied will show a fixed blue light. The timing will re-start to calculate from zero when adjusted during operating.



Intensity adjustment:

Twist the Pressure knob to adjust the pressure intensity, the air pressure intensity range is 50~199mmhg. Rotate to the left, pressure decreases, the lowest air pressure is 50mmhg, rotate to the right, pressure increases, the highest air pressure is 199mmhg (Increase or decrease by 10mmhg). It is recommended to choose the lowest in the first use, then adjust it to personal preference in subsequent use.

Adjusting pressure intensity during treatment does not result in re timing.

Finishing Treatment

The device will automatically shut off after the treatment time is over, or you can press the “Standby” to stop treatment at any time.

	When you end using, pull out the power plug of the main unit from the socket. Caution: when you pull off the power plug or connector, pull it off by holding the plug, not the cord.
	Remove the connector off the main unit
	Remove two attachment connectors from the joint box.
	Keep the main body and all parts in the box.

Device Cleaning Instructions

MAKE SURE TO CUT OFF THE POWER WHEN YOU CLEAN THE MACHINE.
After treatment, wipe the inside and outside of the leg cuffs with clean cotton soft cloth dipped by water.

It's possible to apply 75% alcohol for cleaning the cuffs. Let's the cuffs naturally dry out.

Cautions:

- Don't use gasoline, diluent and other irritating liquid to wipe the machine in case it causes malfunction or the components are damaged or discolored.
- Do not subject the device to moisture or dampness. Do not hold the device under running water, do not submerge it in water or other liquids.
- Ensure that no water penetrates into the machine. If this should ever happen, only use the device again once it is completely dry.
- Not to clean the device during treatment, be sure that the device is turned off when cleaning the device.

Storage

Store the machine in a cool, well-ventilated place. If you wish to keep the device for a long time, place it in its original box.

Do not keep the device at too low temp, or in humid places due to the possibility of freezing and functional damage. Please note the storage conditions as specified in this manual.

Do not store the garments near a stove, cigarette or other heat sources due to the risk of a fire. Do not store the garments near a needles, scissors or other sharp objects due to the risk of damage to the garments.

Avoid the direct sunlight and keep the device out of the heating equipment including gas stove due to the possibility of discoloration and malfunction.
Don't place heavy stuffs on it.

Disposal

At the end of the product life cycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment. Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment.

Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health, you have obligation to dispose of device correctly. Dispose of at public collection point in the EU countries — 2002/96/EC WEEE Directive. Consult your municipal authority or your dealer for information about disposal.

Troubleshoot

Problems	Causes & Solutions
1. Doesn't work and indicator light is off.	1. Check if the device is plugged in. 2. Turn on the power switch. 3. Check if the power cord is malfunctioning, (check up the voltage) 4. Check if the treatment timing is over. (If it is, just press the power button to re-start.)
2. Cuffs do not inflate air or have no pressure.	1. It works only when 2 air hoses are connected properly and completely. 2. Check if hose connectors are twisted or pressed down by heavy objects. 3. Check if all connections of connectors or the power plug are correct and stable.
3. If it sounds like a air leakage.	1. Check if hose connectors are plugged in completely and properly. 2. Check if hose connectors are pressed, twisted or bent. 3. Check if there is any damage to hose connectors
4. Sudden interruption of operating.	1. The power cord or air hoses fall off; 2. The device shuts off automatically after finishing the treatment
5. Pressure is too weak or too strong.	1. Adjust the pressure knob to choose a proper intensity level. 2. Turn off the device if you cannot bear the strength.
6. The main unit gets hot	It is normal if the main unit gets hot after too long time use. We suggest you shut it down for 10 minutes.

Packing list

- 2 x Leg cuffs
- 1 x Joint air hose
- 1 x Main unit
- 1 x User manual

Symbol description

	General warning sign		The product should be vertically up
	Applied part of type BF		Pile Limit 5 layers
	Manufacture date		Standby Button
	Serial Number		Type of protection against electric shock: Class II
	Keep Dry		Manufacturer information: The manufacturer name and address
	Refer to instruction manual.		Batch code
	Hang dry		Do not use iron
	Dry cleaning is not allowed.		Do not bleach
	Upper temperature limits(55°C) and lower limited (-20°C)		
	The range of atmospheric pressure (upper limited : 75kPa , Lower limited : 106kPa) to which the medical device can be safely exposed.		
	The range of humidity (upper limited : 90%RH , Lower limited :5%RH) to which the medical device can be safely exposed.		
	Disposal in accordance with Directive 2002/96/EC (WEEE)		
	This device has not been tested for use in an MR environment and should not be used exposed to MR environments while patients are wearing the device. Keep it outside the MRI scanner room.		

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